

SanipullBulk Premium™ Product Codes

Bulk Sanipull® ribbon for over bed light pulls

Bulk Sanipull cleanable ribbon: spool or coil, standard or anti-ligature¹ option. Each over bed light pull requires a pair of SuperSanipull spring clips²

Order Code	Product Description
SRS1000S-B/G	1000' spooled standard Sanipull ribbon. Blue or Green are standard colours
SRS1000PER-B/G	1000' spooled anti-ligature Sanipull ribbon. Blue or Green are standard colours
SRS360S-B/G	360' spooled standard Sanipull ribbon. Blue or Green are standard colours
SRS360PER-B/G	360' spooled anti-ligature Sanipull ribbon. Blue or Green are standard colours
SRC120S-B/G	120' coiled standard Sanipull ribbon. Blue or Green are standard colours
SRC120PER-B/G	120' coiled anti-ligature Sanipull ribbon. Blue or Green are standard colours



SuperSanipull™ spring clips

Clips are needed on both ends of over bed light pull ribbon. Upper clip attaches to pull chain and lower clip looped around safety rail or onto linen

Order Code	Product Description
SSC10 (10 per bag)	SuperSanipull clip for attachment onto pull chain, linen, or looped around bed rail as shown. White only.



Bulk Sanipull® ribbon for bed side and washroom / tubroom nurse call pulls

Cleanable red ribbon: spool or coil, standard or anti-ligature option. Also available in translucent red (TR), clear (C), orange (O), white (W) or black (K)³


Order Code	Product Description
SRS1000S-R	1000' spooled standard Sanipull red ribbon.
SRS1000PER-R	1000' spooled anti-ligature Sanipull red ribbon.
SRS360S-R	360' spooled standard Sanipull red ribbon.
SRS360PER-R	360' spooled anti-ligature Sanipull red ribbon.
SRC120S-R	120' coiled standard Sanipull red ribbon.
SRC120PER-R	120' coiled anti-ligature Sanipull red ribbon.



SuperSanipull component fittings used with nurse call pulls

Clip attaches to user end of bed side nurse call pull. Adapter adapts to call fixture's string, cord or chain. Handles can be fitted anywhere on the ribbon. Adapter and grab handles also available in white (W) or orange (O).

Order Code	Product Description
SSC10 (10 per bag)	Spring clip for attachment to linen. White clips only.
SSNCEA5 (5 per bag)	Modified clamshell adapter for making temporary bed side nurse call pull extension. Red is standard
SSA10-R (10 per bag)	Clamshell adapter for fast and esthetic installation. Red is standard
SSH10-R (10 per bag)	Grab handle for accessibility anywhere on ribbon. Red is standard



Tools to Assist Installation

Hand punch is needed for standard ribbon and helpful for anti-ligature ribbon installation. Set offered with scissors for convenience.

Order Code	Product Description
PSS	bandage scissors to cut ribbon and medium duty 1/4" hand punch to make hardware attachment holes



Important Product Notes

SanipullBulk is best suited for facilities or departments with more than 140 beds. For complete installation instructions visit www.Sanipull.com

¹The Sanipull anti-ligature ribbon (symbol, PER) option helps mitigate the risk of entanglement or strangulation with a proprietary safety breakaway designed to tear apart at 27 +/- 2 lbs force at 22°C. This ribbon is factory modified with 5/16" diameter holes punched every 6 inches along its length. Because of the nature of this ribbon modification, warranty extends only to the attachment hardware. Note: Wilshire Works or its distributors will not be held liable for death or damages directly or indirectly involving any Sanipull product.

²SuperSanipull spring clips are designed to release from an over bed light pull chain at 12 +/- 1 lbs overforce pull. Simply re-attach for ongoing switch protection. If the user-end of the ribbon is looped around bed rail with the clip snapped onto the ribbon, the clip will release if the safety rail is lowered or the bed is moved.

³Standard Sanipull ribbon colours: blue (B) and green (G) for over bed light pulls and red (R), which is recommended for nurse call pulls. Special order colours: black (K), white (W), orange (O), clear (C) and translucent red (R). Add 10¢/ft.



Sanipull is a durable solution that can be repaired by replacing damaged components. Wilshire Works offers a 5-year replacement guarantee on SuperSanipull components if correctly assembled and installed.

Pre-cut and anchor hole punched Sanipull ribbon (on 1000 ft spools only) for standard length pulls is available. A \$25/spool set up fee will apply. Call for details.

Wilshire Works Inc. continually explore ways to improve our products and therefore reserve the right to make design or material changes. For more information visit www.Sanipull.com.