

2014 Sanipull® Replacement Parts End-User Product Codes

All plastic fittings are warranted for 5-years



SuperSanipull™ Spring Clip

(Releases from fixture ball chain at 12 +/- 2 lb of pull force)
 Multi-purpose spring clip for switch saver attachment to light switch pull chain or as a positive-grip, non-destructive attachment to linen or around a bed rail.

Guaranteed to work and not break under normal working conditions.
 Will not tear bed linen.

Bag of 10 clips

PN: SSC10



SuperSanipull™ Accessibility Handle and Cover

Attaches anywhere on Sanipull ribbon with cover to facilitate easy wipe down. Colour Options: red (R), white (W) or orange (O)

Bag of 10 sets (OPEN handle)

PN: SBPHO10

Bag of 10 sets (CLOSED handle)

PN: SBPHC10



SuperSanipull™ Adapter for Nurse Calls

Provides quick retrofit attachment of Sanipull ribbon to string, nylon cord, ball chain, beaded cord or plastic pull cords. Perfect adapter for pull replacement for Austco™ nylon cord and plastic pull cords. Colour options: red (R), white (W) or orange (O).

Bag of 10 adapters

PN: SBPA10



Bandage Scissors

The right tool you need to quickly install SuperSanipull ribbon pulls

Scissors

PN: SCI

Ball Chain Lengths and Connectors

Use to extend the length of pull switch ball chain to obtain optimal length for over bed Sanipull

10' of #6 nickel plated ball chain and 20 connectors

PN: BCL10

100' of #6 nickel plated ball chain and 200 connectors

PN: BCL100

Bag of 100 metal ball chain connectors

PN: BCC100



Sanipull Poster for staff lounge or nursing station

8 1/2" x 11" full-colour on poster card stock reminds staff to attach Sanipull properly and report problems to maintenance ASAP. Helps save maintenance costs. No cost printable pdf also available on request by email.

5 copies Sanipull Original

PN: SOP5

5 copies SBP & SS RTU

PN: SBPP5

