



2020 Sanipull[™] Overbed Light Pull Product Codes

NEW ** Lifetime warranty on plastic fittings

1-1. SuperSanipull RTU : pre-cut, assembled, Ready To Use

SuperSanipull RTU for over bed light pulls

Pre-cut Sanipull cleanable ribbon pull is assembled with one Sanipull release ball chain adapter and 2 durable Universal Sanipull clips¹, coiled, 10 pulls in a bundle. 3 pull lengths available with standard or anti-ligature² modified ribbon option. Blue colour ribbon.

(SuperSanipull RTU is sold in bundles of 10 units. Below is pull unit price.)

Order Code

Product Description

SSLTRB4SSuperSanipull pre-cut 4 ft standard ribbon over bed light pull with one release ball chain adapter and two clipsSSLTRB4PERSuperSanipull pre-cut 4 ft anti-ligatured ribbon over bed light pull with one release ball chain adapter and two clipsSSLTRB5SSuperSanipull pre-cut 5 ft standard ribbon over bed light pull with one release ball chain adapter and two clipsSSLTRB5PERSuperSanipull pre-cut 5 ft anti-ligatured ribbon over bed light pull with one release ball chain adapter and two clipsSSLTRB6SSuperSanipull pre-cut 5 ft anti-ligatured ribbon over bed light pull with one release ball chain adapter and two clipsSSLTRB6SSuperSanipull pre-cut 6 ft standard ribbon over bed light pull with one release ball chain adapter and two clipsSSLTRB6PERSuperSanipull pre-cut 6 ft anti-ligatured ribbon over bed light pull with one release ball chain adapter and two clipsSSLTRB6PERSuperSanipull pre-cut 6 ft anti-ligatured ribbon over bed light pull with one release ball chain adapter and two clips

1-2. SanipullBulk Premium Kit : DIY bulk, cut any length, assemble on site (all fittings included)

SanipullBulk Premium Kit for over bed light pulls

Kit comes with Sanipull blue ribbon, Sanipull release ball chain adapters, durable Universal Sanipull clips. 2 quantity options are available. Blue colour ribbon is standard. ** Pull quality is based on assembling 6 ft length. Please see the required fitting chart for additional pulls. **

Order Code

Product Description

300 ft Kit (50 pulls Kit)

SBPKLT300S 300 ft of Sanipull standard blue ribbon, 50 Sanipull release ball chain adapters and 100 Universal Sanipull clips SBPKLT300PER 300 ft of Sanipull anti-ligature blue ribbon, 50 Sanipull release ball chain adapters and 100 Universal SuperSanipull clips

900 ft Kit (150 pulls Kit)

SBPKLT900S 900 ft of Sanipull standard blue ribbon, 150 Sanipull release ball chain adapters and 300 Universal SuperSanipull clips SBPKLT900PER 900 ft of Sanipull anti-ligature blue ribbon, 150 Sanipull release ball chain adapters and 300 Universal SuperSanipull clips

Tools to Assist Installation

Scissors	Serrated edge bandage scissors to cut ribbon
Sanipull fitting pry blade	Fine profile blade that makes opening a clip or adapter clamshell door a snap

1-3. Sanipull fittings for over bed light pulls

*Fixture end of the pull will require both release ball chain adapter and universal Sanipull clip		F
Order Code	Product Description	N
USC10	Bag of 10 Universal Sanipull spring clips. This multipurpose clip used on over bed lights is attached onto ball chain adapter (below) to protect fixture from over force pull damage. Also, clip on lower end of pull allows convenient positioning.	Dates
BCA10	Bag of 10 Sanipull release ball adapters. Securely connects to fixture switch ball chain. Sanipull clip attaches to the ball end.	

Important Product Notes

For complete installation instructions or support information visit www.Sanipull.com or call 1-800-SANIPUL

¹Sanipull spring clip attached onto ball chain adapter provides consistent and reliable 13 +/- 1 lbs (5 +/- 0.5 kg) overforce pull protection for over bed light fixtures. Simply re-attach for ongoing switch protection. If the user-end of the ribbon is attached onto linen or looped around bed rail with the clip snapped onto the ribbon, the clip will harmlessly release if the safety rail is lowered or the bed is moved.

²Sanipull anti-ligature ribbon (symbol, PER) option helps mitigate the risk of entanglement or strangulation with a proprietary safety breakaway designed to tear apart at 29 +/- 2 lbs (13 +/-1 kg) force at 74 °F (23 °C). Anti-ligature ribbon is factory modified with nominal 5/16" diameter holes punched every 6-7 inches along its length. Because of the nature of this ribbon modification, warranty extends only to the attachment fittings. **Note:** Wilshire Works or its distributors will not be held liable for death or damages directly or indirectly involving any Sanipull product.

Standard Sanipull ribbon colours: blue (B) for over bed light pulls and red (R) for nurse call pulls.

Sanipull is a durable solution that can be repaired by replacing damaged components. Wilshire Works offers a 100% end user satisfaction guarantee and a lifetime replacement warranty on SuperSanipull fittings if correctly installed.

Wilshire Works continually works to improve its products and therefore reserves the right to make design or material changes.





